


## 小分子-靶点蛋白富集及质谱筛选 (**LFQ**)

SM-PD Enrichment Kit (**Biotin** tag, MS grade)

Catalog Number: MG05

## Workflow of SM-PD Enrichment Kit (Biotin tag, MS grade) Catalog Number: #MG05

### PART 1 Guide for large-scale cell culture, protein extraction & Preclear



1. Collect cell pellets
2. Liquid Nitrogen snap-freezing
3. Store at -80°C

1. Extract proteins with **MG05-1**
2. Filter protein solution
3. Protein quantitation (Imultiomics, #MGR02)

Protein preclear with **MG05-5**

Store at -80°C or go to downstream applications

### PART 2 **Probe** concentration optimization

Start at least **21 mg**  
Protein sample  
(**MG05-1** extracted)

1. Preclear with **MG05-5**
2. Divide into **5** parts

Perform SM-binding reactions with **MG05-4** and **Probe**

Perform pull-down with **MG05-5**

Wash with **MG05-2 & MG05-3**

Elute with **MG05-6**

SDS-PAGE & staining (Imultiomics, #MGR01)

1. Obtain optimized **Probe** concentration
2. Obtain optimized **Protein input**/reaction (**X mg**)
3. Obtain optimized **MG05-5** input/reaction (**μl**)

\***Probe**: Biotin-labeled SM (**User Provide**)

### PART 3 **SM** concentration optimization

at least **6× Protein input (X mg)**  
Protein sample  
(**MG05-1** extracted)

1. Preclear with **MG05-5**
2. Divide into **6** parts

Perform SM-binding reactions with **MG05-4**, **SM**, and **optimized Probe**

Perform pull-down with optimized volume (**μl**) of **MG05-5**

Wash with **MG05-2 & MG05-3**

Elute with **MG05-6**

SDS-PAGE & staining (Imultiomics, #MGR01)

Obtain optimized **SM** concentration

\***SM**: parent compound (**User Provide**)

### PART 4 SM-PD sample preparation for MS detection

at least **3× Protein input (X mg)**  
Protein sample  
(**MG05-1** extracted)

1. Preclear with **MG05-5**
2. Divide into **3** parts

Perform SM-binding reactions with **optimized Probe & SM**

Perform pull-down with optimized volume (**μl**) of **MG05-5**


Wash with **MG05-2 & MG05-3**

On-beads trypsin digestion  
(Imultiomics, #MG04)

LC-MS detection & LFQ quantification

Statistical analysis to distinguish **SM-targeted Proteins**

A. Protein extraction, Quantitation & Preclear




1. Cell pellet of 7×10cm dishes were lysed with 7×600μl of **MG05-1** (Cell lysis buffer).
2. Protein concentration (**6.8 mg/ml**) was quantified using Bradford Protein Assay Kit (**Imultiomics, #MGR02**).
3. Protein sample was precleared with MG05-5 according to the protocol of MG05 Kit

曲线名称	曲线公式	A	B	R2
标准曲线	Y=A*X+B	0.0738	1.2	0.991




## B. Pull-down data of PC in MG05 Kit



1. BP: Before Preclear
2. AP: After Preclear
3. NC: negative control (MG05-5 beads only)
4. PC (Positive control): MG05-4 in **MG05 Kit**. PC is the biotin-labeled small molecule (SM).  
The pull-down data of PC showed differential bands in contrast to that of NC.
5. SDS-PAGE gel was stained with Coomassie Blue Staining Reagent (**Imultiomics, #MGR01**)

### C. Streptavidin remained on beads after on-beads digestion



**PD:** MG05-5 (Streptavidin-beads) enriched proteins by pull-down. PD was conducted in duplicate.

**A.** PD sample was carried out on-beads digestion using Co-IP/pull-down in-solution digestion Kit (**Imultiomics, #MG04**), the enriched proteins were completely digested into peptides in contrast to that of PD.

**B.** The left beads after on-beads digestion was treated with MG05-6 (5× sample loading buffer) for Coomassie blue staining (**Imultiomics, #MGR01**). The data showed most of streptavidin remained on beads (MG05-5) after on-beads digestion, less interference of Pull-down enriched proteins identification by LC-MS!